

Anti-SCAMP5 antibody (180-230 aa) (STJ13102030)

STJ13102030

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SCAMP5 (180-230 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	56807
Gene Symbol	Scamp5
Uniprot ID	SCAM5_MOUSE
Immunogen	A synthetic peptide from aa region 180-230 of rat SCAMP5 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in mouse and shares 95% identity with human NA s.
Immunogen Region	180-230 aa
Specificity	Specific for SCAMP5.
Immunogen Sequence	