

Anti-SLC16A2 antibody (10-60 aa) (STJ13102151)

STJ13102151

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC16A2 (10-60 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	6567
Gene Symbol	SLC16A2
Uniprot ID	MOT8_HUMAN
Immunogen	A synthetic peptide from aa region 10-60 of human SLC16A2 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	10-60 aa
Specificity	Specific for SLC16A2.
Immunogen Sequence	