

Anti-SLC3A2 antibody (Internal Region) (STJ13102165)

STJ13102165

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC3A2 (Internal Region) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	6520
Gene Symbol	SLC3A2
Uniprot ID	4F2_HUMAN
Immunogen	A synthetic peptide from internal region of human CD98 (SLC3A2, 4F2 heavy chain antigen, MDU1) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	Internal Region
Specificity	Specific for SLC3A2.
Immunogen Sequence	