

Anti-SLC45A1 antibody (C-Term) (STJ13102179)

STJ13102179

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC45A1 (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). A dilution of 1: 200 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	50651
Gene Symbol	SLC45A1
Uniprot ID	S45A1_HUMAN
Immunogen	A synthetic peptide from c-terminal region of human SLC45A1 (Proton-associated sugar transporter A, PAST-A) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	C-Term
Specificity	Specific for SLC45A1.
Immunogen Sequence	