

Anti-SLC30A7 antibody (200-250 aa) (STJ13102226)

STJ13102226

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-SLC30A7 (200-250 aa) is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC/WB |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | IgG purified |
| Dilution Range | IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications. |
| Formulation | Lyophilised |
| Isotype | IgG |
| Storage | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. |
| Instruction | When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 148867 |
| Gene Symbol | SLC30A7 |
| Uniprot ID | ZNT7_HUMAN |
| Immunogen | A synthetic peptide from aa region 200-250 of human SLC30A7 conjugated to an immunogenic carrier protein was used as the antigen. |
| Immunogen Region | 200-250 aa |
| Specificity | Specific for SLC30A7. |
| Immunogen Sequence | |