

**Anti-SLC9A8 antibody (450-500 aa) (STJ13102278)**

STJ13102278

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-SLC9A8 (450-500 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	IgG purified
<b>Dilution Range</b>	IHC, WB (confirmed by recombinant protein). A dilution of 1: 200 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
<b>Formulation</b>	Lyophilised
<b>Isotype</b>	IgG
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
<b>Instruction</b>	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">23315</a>
<b>Gene Symbol</b>	<a href="#">SLC9A8</a>
<b>Uniprot ID</b>	<a href="#">SL9A8_HUMAN</a>
<b>Immunogen</b>	A synthetic peptide from aa region of 450-500 of human SLC9A8 conjugated to blue carrier protein was used as the antigen. The peptide shares 88% identity with mouse and rat sequences.
<b>Immunogen Region</b>	450-500 aa
<b>Specificity</b>	Specific for SLC9A8.
<b>Immunogen Sequence</b>	