

**Anti-SLC17A9 antibody (200-250 aa) (STJ13102282)**

STJ13102282

**GENERAL INFORMATION**

|                          |                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                       |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-SLC17A9 (200-250 aa) is suitable for use in Immunohistochemistry and Western Blot research applications. |
| <b>Applications</b>      | IHC/WB                                                                                                                                   |
| <b>Host/Source</b>       | Rabbit                                                                                                                                   |
| <b>Reactivity</b>        | Human                                                                                                                                    |

**PRODUCT PROPERTIES**

|                       |                                                                                                                                                                                              |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Polyclonal                                                                                                                                                                                   |
| <b>Clone ID</b>       |                                                                                                                                                                                              |
| <b>Concentration</b>  |                                                                                                                                                                                              |
| <b>Conjugation</b>    | Unconjugated                                                                                                                                                                                 |
| <b>Purification</b>   | IgG purified                                                                                                                                                                                 |
| <b>Dilution Range</b> | IHC, WB (confirmed by recombinant protein). A dilution of 1: 200 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications. |
| <b>Formulation</b>    | Lyophilised                                                                                                                                                                                  |
| <b>Isotype</b>        | IgG                                                                                                                                                                                          |
| <b>Storage</b>        | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.                                                        |
| <b>Instruction</b>    | When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.                                                                                  |

**TARGET INFORMATION**

|                           |                                                                                                                         |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">63910</a>                                                                                                   |
| <b>Gene Symbol</b>        | <a href="#">SLC17A9</a>                                                                                                 |
| <b>Uniprot ID</b>         | <a href="#">S17A9_HUMAN</a>                                                                                             |
| <b>Immunogen</b>          | A synthetic peptide from aa region 200-250 of human SLC17A9 conjugated to blue carrier protein was used as the antigen. |
| <b>Immunogen Region</b>   | 200-250 aa                                                                                                              |
| <b>Specificity</b>        | Specific for SLC17A9.                                                                                                   |
| <b>Immunogen Sequence</b> |                                                                                                                         |