

Anti-TRPC2 antibody (STJ13102312)

STJ13102312

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-TRPC2 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Sheep
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. A dilution of 1: 1000 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	22064
Gene Symbol	Trpc2
Uniprot ID	TRPC2_MOUSE
Immunogen	A synthetic peptide from mouse TRPC2 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in rat.
Region	
Specificity	Specific for TRPC2.
Immunogen Sequence	