

Anti-TrkB antibody (400-450 aa) (STJ13102380)

STJ13102380

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-TrkB (400-450 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Sheep
Reactivity	Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A concentration of 10-50 ug/ml is recommended. The optimal concentration should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A synthetic peptide from aa region 400-450 of the extracellular domain of mouse TrkB conjugated to blue carrier protein was used as the antigen.
Immunogen Region	400-450 aa
Specificity	Specific for TrkB.
Immunogen Sequence	