

Anti-Phospho-TrkC-Y820 antibody (STJ13102446)

STJ13102446

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-TrkC-Y820 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Rat/Mouse/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:3000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications. This antibody may also recognise the unphosphorylated form.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A synthetic peptide to phospho Y820 of human TrkC conjugated to blue carrier protein was used as the antigen. The antigen is homologous in rat and mouse.
Immunogen Region	
Specificity	Specific for TrkC.
Immunogen Sequence	