

Anti-TSG101 antibody (280-330 aa) (STJ13102465)

STJ13102465

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TSG101 (280-330 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	22088
Gene Symbol	Tsg101
Uniprot ID	TS101_MOUSE
Immunogen	A synthetic peptide from aa region 280-330 of rat TSG101 conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	280-330 aa
Specificity	Specific for TSG101.
Immunogen Sequence	