

Anti-TAS2R31 antibody (260-309 aa) (STJ13102565)

STJ13102565

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAS2R31 (260-309 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	259290
Gene Symbol	TAS2R31
Uniprot ID	T2R31_HUMAN
Immunogen	A synthetic peptide from aa region 260-309 of human TAS2R31 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	260-309 aa
Specificity	Specific for TAS2R31 but some cross-reactivity with T2R46 and T2R43 is expected.
Immunogen Sequence	