

Anti-TAS2R50 antibody (200-250 aa) (STJ13102569)

STJ13102569

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAS2R50 (200-250 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	259296
Gene Symbol	TAS2R50
Uniprot ID	T2R50_HUMAN
Immunogen	A synthetic peptide from aa region 200-250 of human TAS2R50 conjugated to blue carrier protein was used as the antigen. The peptide shares about 90% with TAS2R43, TAS2R31, TAS2R30, TAS2R20.
Immunogen Region	200-250 aa
Specificity	Specific for TAS2R50.
Immunogen Sequence	