

Anti-TAS2R7 antibody (260-318 aa) (STJ13102576)

STJ13102576

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAS2R7 (260-318 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	50837
Gene Symbol	TAS2R7
Uniprot ID	TA2R7_HUMAN
Immunogen	A synthetic peptide from aa region 260-318 of human TAS2R7 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	260-318 aa
Specificity	Specific for TAS2R7.
Immunogen Sequence	