

Anti-ULK4 antibody (STJ13102593)

STJ13102593

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ULK4 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). Use at a dilution of 1:1000. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	54986
Gene Symbol	ULK4
Uniprot ID	ULK4_HUMAN
Immunogen	A synthetic peptide from human ULK4 conjugated to blue carrier protein has been used as the antigen. The peptide shares 95% identity with mouse sequence.
Immunogen Region	
Specificity	Specific for ULK4.
Immunogen Sequence	