

## Anti-SNF8 antibody (Internal Region) (STJ13102608)

STJ13102608

### GENERAL INFORMATION

|                          |                                                                                                                                            |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                         |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-SNF8 (Internal Region) is suitable for use in Immunohistochemistry and Western Blot research applications. |
| <b>Applications</b>      | IHC/WB                                                                                                                                     |
| <b>Host/Source</b>       | Rabbit                                                                                                                                     |
| <b>Reactivity</b>        | Mouse/Rat/Human                                                                                                                            |

### PRODUCT PROPERTIES

|                       |                                                                                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Polyclonal                                                                                                                                              |
| <b>Clone ID</b>       |                                                                                                                                                         |
| <b>Concentration</b>  |                                                                                                                                                         |
| <b>Conjugation</b>    | Unconjugated                                                                                                                                            |
| <b>Purification</b>   | IgG purified                                                                                                                                            |
| <b>Dilution Range</b> | IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications. |
| <b>Formulation</b>    | Lyophilised                                                                                                                                             |
| <b>Isotype</b>        | IgG                                                                                                                                                     |
| <b>Storage</b>        | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.                   |
| <b>Instruction</b>    | When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.                                             |

### TARGET INFORMATION

|                           |                                                                                                                                                                                       |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 27681                                                                                                                                                                                 |
| <b>Gene Symbol</b>        | Snf8                                                                                                                                                                                  |
| <b>Uniprot ID</b>         | SNF8_MOUSE                                                                                                                                                                            |
| <b>Immunogen</b>          | A synthetic peptide from the internal region of human SNF8 conjugated to an immunogenic carrier protein was used as the antigen. The antigen is homologous in rat, mouse and xenopus. |
| <b>Immunogen Region</b>   | Internal Region                                                                                                                                                                       |
| <b>Specificity</b>        | Specific for SNF8.                                                                                                                                                                    |
| <b>Immunogen Sequence</b> |                                                                                                                                                                                       |