

Anti-VAMP4 antibody (STJ13102638)

STJ13102638

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-VAMP4 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Goat
Reactivity	Human/Mouse/Rat/Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	8674
Gene Symbol	VAMP4
Uniprot ID	VAMP4_HUMAN
Immunogen	A synthetic peptide from mouse VAMP4 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for VAMP4.
Immunogen Sequence	