

Anti-GST antibody (STJ140007)

STJ140007

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Gluthathione-S-Transferase is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	GST

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000 Catarino S Ribeiro-Rodrigues TM Sa Ferreira R et al. Cells 2020 Apr. PMID: 32272685 Seixas C Choi SY Polgar N et al. Mol Biol Cell. 2015 Nov 18. PMID: 26582389
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified recombinant protein produced in E. coli.
Immunogen Region	
Specificity	In lysates of E. coli cells transformed with pGEX plasmid detects a band of 26 kDa by Western blot.
Immunogen Sequence	



Anti-GST antibody at 1:1000 dilution; 50 ng of protein per lane; rabbit polyclonal to goat IgG (HRP) at 1:10000 dilution

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081