

Anti-RYMV P1 antibody (55aa to N-terminal) (STJ140101)

STJ140101

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-RYMV P1 (55aa to N-terminal) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	Transfected cells

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified recombinant peptide derived from within residues 55 aa to N-terminal of RYMV P1 produced in E. coli.
Immunogen Region	55aa to N-terminal
Specificity	Detects RYMV P1 full length recombinant fusion protein GFP-P1 expressed in transfected mammalian cells.
Immunogen Sequence	



Anti-RYMV P1 antibody at 1:1000 dilution; MBP-RYMV P1 N-term recombinant protein and HEK293 transfected cell lysate; rabbit polyclonal to goat IgG (HRP) at 1:10000 dilution.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081