

Anti-SNCG antibody (80aa C-Term) (STJ140141)

STJ140141

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Synuclein gamma (80aa C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	Human/Rat/Mouse/Monkey/Canine

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:5000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	6623
Gene Symbol	SNCG
Uniprot ID	SYUG_HUMAN
Immunogen	Purified recombinant peptide derived from within residues 80 aa to the C-terminus of human SNCG produced in E. coli.
Immunogen Region	80aa C-Term
Specificity	Detects endogenous levels of SNCG by Western blot. This Ab is specific for SNCG and does not recognizes SNCA and SNCB proteins.
Immunogen Sequence	



Anti-SNCG antibody at 1:2500 dilution; lysate at 50 Åug per lane; rabbit polyclonal to goat IgG (HRP) at 1:10000 dilution