

Anti-CHML antibody (STJ140146)

STJ140146

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Choroideremia (Rab escort protein 2) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	Human/Rat/Mouse/Monkey/Canine

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:5000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	1122
Gene Symbol	CHML
Uniprot ID	RAE2_HUMAN
Immunogen	Purified recombinant human TNF alpha produced in E. coli.
Immunogen Region	
Specificity	Detects endogenous REP2 by Western blot. This Ab is specific for REP2, it does not cross-react with REP1.
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081