

Anti-TBX18 antibody (330aa C-Term) (STJ140163)

STJ140163

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-T-box transcription factor TBX18 (330aa C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	Human/Rat/Mouse/Monkey/Canine

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	9096
Gene Symbol	TBX18
Uniprot ID	TBX18_HUMAN
Immunogen	Recombinant peptide derived from within residues 330 aa to the C-terminus of Tbx18 produced in E. coli.
Immunogen Region	330aa C-Term
Specificity	Detects endogenous levels of TBX18 in AtT-20, MCF7, HeLa, H69 and Brain by Western blot.
Immunogen Sequence	



Anti-TBX18 antibody at 1:1000 dilution; lysates at 50 µg per lane; rabbit polyclonal to goat IgG (HRP) at 1:10000 dilution

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081