

Anti-Streptavidin antibody (50 aa C-Term) (STJ140171)

STJ140171

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Streptavidin (50 aa C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	RecombinantStreptavidin

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:5000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Peptide derived from within residues 50 aa to the C-terminus of streptavidin protein.
Immunogen Region	50 aa C-Term
Specificity	Reacts specifically with streptavidin protein.
Immunogen Sequence	CPSAASIDAACKAGVNNGNP LDAVQQ



Anti-Streptavidin antibody at 1/2, 500 dilution; lysates contain 30 ng of recombinant protein; rabbit polyclonal to goat IgG (HRP) at 1:10000 dilution

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081