

Anti-SARS-CoV-2 Membrane Protein antibody (100-stop) (STJ140196)

STJ140196

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-SARS-CoV-2 Membrane Protein (100-stop) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Affinity purified recombinant fusion protein using the C-terminal of Membrane protein (residues 100 to stop) and produced in E. coli.
Immunogen Region	100-stop
Specificity	In lysates of transfected cells with the plasmid containing the sequence used, detects the fusion protein by Western blot.
Immunogen Sequence	



Anti-Membrane protein antibody at 1/2, 500 dilution;
lane with 30 ng of recombinant fusion protein (103-stop
aa) ; rabbit polyclonal to goat IgG (HRP) at 1:10000
dilution

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081