

Anti-SARS-CoV-2 Spike-RBD Domain antibody (343-400) (STJ140197)

STJ140197

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-SARS-CoV-2 Spike-RBD Domain (343-400) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Affinity purified recombinant fusion protein using the spike protein RBD domain (residues 343 to 400) and produced in E. coli.
Immunogen Region	343-400
Specificity	In lysates of transfected cells with the plasmid containing the sequence used, detects the fusion protein by Western blot.
Immunogen Sequence	



Anti-Spike antibody at 1/2, 500 dilution; lane with 30 ng of recombinant fusion protein; chicken polyclonal to goat IgG (HRP) at 1:10000 dilution

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081