

Goat Anti-Rabbit IgG (Full length) antibody {DyLight®405} (STJ140238)

STJ140238

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	DyLight®405-conjugated goat polyclonal anti-Rabbit IgG Full length secondary antibody. For use in most research applications.
Applications	WB/IF
Host/Source	Goat
Reactivity	Rabbit

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	DyLight®405
Purification	Epitope affinity purified
Dilution Range	WB 1:2000-1:10000 IF 1:100-1:1000
Formulation	PBS, 20% Glycerol, 0.05% ProClin® and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Full rabbit IgG affinity purified from domestic rabbit (O. coniculus) serum using protein A. Purity >95% by SDS-PAGE.
Immunogen Region	Full length
Specificity	Using serum samples detects rabbit IgG by Western blot. Reduced cross-reactivity to other species.
Immunogen Sequence	