

Goat Anti-Human IgA (Full length) antibody {Biotin} (STJ140249)

STJ140249

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated goat polyclonal anti-Human IgA Full length secondary antibody. For use in most research applications.
Applications	WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Biotin
Purification	Epitope affinity purified
Dilution Range	WB 1:2000-1:10000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Full human IgA affinity purified (purity >95% by SDS-PAGE).
Immunogen Region	Full length
Specificity	Using serum samples detects the light and heavy chains by Western blot.
Immunogen Sequence	