

Anti-GAPDH antibody (240 aa C-Term) {DyLight®633} (STJ140267)

STJ140267

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-GAPDH (240 aa C-Term) is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC-F/IF
Host/Source	Goat
Reactivity	Human/Rat/Mouse/Zebrafish/Canine/Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	DyLight®633
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:5000
	IHC-F 1:200-1:1000
	IF 1:50-1:250
Formulation	PBS, 20% Glycerol, 0.05% ProClin® and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	2597
Gene Symbol	GAPDH
Uniprot ID	G3P_HUMAN
Immunogen	Purified recombinant peptide derived from within residues 240 aa to the C-terminus of human GAPDH produced in E. coli.
Immunogen Region	240 aa C-Term
Specificity	Detects a band of 37 kDa by Western blot.
Immunogen Sequence	