

Anti-His-Tag antibody (STJ140291)

STJ140291

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-His-Tag is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	Reacts against HIS tagged proteins

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB:1:1, 000-1:5, 000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified recombinant peptide produced in E. coli.
Immunogen Region	
Specificity	The antibody recognizes the His-tag fused to the amino-or carboxy-termini of targeted proteins in transformed or transfected cells.
Immunogen Sequence	