

**Anti-GAPDH antibody (C-Term) {DyLight®755} (STJ140311)**

STJ140311

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Goat polyclonal antibody anti-GAPDH (C-Term) is suitable for use in Western Blot, Immunofluorescence and Immunohistochemistry research applications.
<b>Applications</b>	WB/IF/IHC-F
<b>Host/Source</b>	Goat
<b>Reactivity</b>	Reacts against human, rat, mouse, zebrafish, canine and monkey proteins.

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	2.5 mg/mL
<b>Conjugation</b>	DyLight®755
<b>Purification</b>	Epitope affinity purified
<b>Dilution Range</b>	WB:1:500-1:5,000
	IF:1:50-1:250
	IHC-F:1:200-1:1,000
<b>Formulation</b>	PBS, 20% Glycerol and 0.05% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Purified recombinant peptide derived from within residues 240 aa to the C-terminus of human GAPDH produced in E. coli.
<b>Immunogen Region</b>	C-Term
<b>Specificity</b>	Detects a band of 37 kDa by Western blot.
<b>Immunogen Sequence</b>	