

Anti-CD3 antibody [SPV-T3b] (STJ16100004)

STJ16100004

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD3 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-F/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SPV-T3b
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence