

Anti-CD8 antibody [FK18] (STJ16100010)

STJ16100010

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD8 is suitable for use in Immunofluorescence research applications.
Applications	IF
Host/Source	Mouse
Reactivity	Human/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FK18
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IF 1:10
Formulation	Serum free medium containing 1% FCS
Isotype	IgG3
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human thymocytes
Immunogen Region	
Specificity	
Immunogen Sequence	