

Anti-CD57 antibody [6-13-19-1] (STJ16100011)

STJ16100011

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD57 is suitable for use in Immunofluorescence and Immunohistochemistry research applications.
Applications	IF/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6-13-19-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IF 1:5 IHC 1:5
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	JURKAT T-cell line
Immunogen Region	
Specificity	
Immunogen Sequence	