

Anti-HLA-DRw52 antibody [7.3.19.1] (STJ16100015)

STJ16100015

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-DRw52 is suitable for use in Immunofluorescence research applications.
Applications	IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	7.3.19.1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IF 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human HLA-DRw52
Immunogen Region	
Specificity	
Immunogen Sequence	