

Anti-CD26 antibody [202.36] (STJ16100026)

STJ16100026

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD26 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	FC/IHC-F/IF/WB
Host/Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	202.36
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 0 5-1 0 µg/million cells in 0 1 mL-WB 1-2 µg/mL-IF 0 5-1 0 µg/mL-IHC Re
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1803
Gene Symbol	DPP4
Uniprot ID	DPP4_HUMAN
Immunogen	T-cell clone
Immunogen Region	
Specificity	
Immunogen Sequence	