

Anti-CD37 antibody [IPO-24] (STJ16100030)

STJ16100030

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD37 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC/IHC-F/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IPO-24
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 0 5-1 0 µg/million cells in 0 1 mL-IF 0 5-1 0 µg/mL-IHC 1-2 µg/ml for 30 min at RT
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	951
Gene Symbol	CD37
Uniprot ID	CD37_HUMAN
Immunogen	spleen cells of a patient with hairy cell leukemia
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	