

Anti-CD22 antibody [RFB4] {FITC} (STJ16100031)

STJ16100031

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD22 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RFB4
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 4mg/ml BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	933
Gene Symbol	CD22
Uniprot ID	CD22_HUMAN
Immunogen	Human CD22
Immunogen Region	
Specificity	
Immunogen Sequence	