

Anti-CD147 antibody [MEM-M6/1] (STJ16100037)

STJ16100037

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD147 is suitable for use in Flow Cytometry, Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	FC/IP/WB/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-M6/1
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Protein A-CR purified soluble recombinant form of CD147, CD147Rg, which consists of the cDNA coding for the hinge region, CH2- and CH3 domain of human IgG1 (CD147Rg is secreted by transfectants as a dimer).
Immunogen Region	
Specificity	
Immunogen Sequence	