

Anti-HLA-DQ1+DQ3 antibody [HL-37] (STJ16100048)

STJ16100048

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-DQ1+DQ3 is suitable for use in Flow Cytometry, Immunoprecipitation and Western Blot research applications.
Applications	FC/IP/WB
Host/Source	Mouse
Reactivity	Pig/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HL-37
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 5 µg/ml-IP-WB to be determined
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG3
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Burkitt's lymphoma cell line Raji.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081