

Anti-CD8 antibody [MEM-31] (STJ16100051)

STJ16100051

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD8 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-31
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1 µg/ml-IP-Mass Cytometry
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Crude thymus membrane fraction.
Immunogen Region	
Specificity	
Immunogen Sequence	