

Anti-CALLA antibody [MEM-78] (STJ16100053)

STJ16100053

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD10/CALLA is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-78
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1 µg/ml-IP
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4311
Gene Symbol	MME
Uniprot ID	NEP_HUMAN
Immunogen	NALM-6 human pre-B cell line
Immunogen Region	
Specificity	
Immunogen Sequence	