

Anti-CD11a antibody [MEM-25] (STJ16100054)

STJ16100054

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD11a is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-25
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 2 µg/ml-IP excellent antibody for immunoaffinity purification of LFA-1 complex
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3683
Gene Symbol	ITGAL
Uniprot ID	ITAL_HUMAN
Immunogen	Leukocytes from a patient suffering from a LGL-type leukaemia.
Immunogen Region	
Specificity	
Immunogen Sequence	