

Anti-CD20 antibody [MEM-97] (STJ16100060)

STJ16100060

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD20 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
Applications	FC/IP
Host/Source	Mouse
Reactivity	Pig/Cow/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-97
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 10 µg/ml-IP
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	931
Gene Symbol	MS4A1
Uniprot ID	CD20_HUMAN
Immunogen	Raji human Burkitt's lymphoma cell line
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081