

Anti-CD53 antibody [MEM-53] (STJ16100066)

STJ16100066

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD53 is suitable for use in Flow Cytometry, Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	FC/IP/WB/IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-53
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	963
Gene Symbol	CD53
Uniprot ID	CD53_HUMAN
Immunogen	Leukocytes of patient suffering from a LGL-type leukemia.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	