

Anti-CD14 antibody [Tük 4] (STJ16100069)
STJ16100069**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD14 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Non-Human Primates/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Tük 4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	929
Gene Symbol	CD14
Uniprot ID	CD14_HUMAN
Immunogen	Human CD14
Immunogen Region	
Specificity	
Immunogen Sequence	