

Anti-CD59 antibody [MEM-43] (STJ16100070)

STJ16100070

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD59 is suitable for use in Flow Cytometry, Immunoprecipitation and Immunohistochemistry research applications.
Applications	FC/IP/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-43
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	966
Gene Symbol	CD59
Uniprot ID	CD59_HUMAN
Immunogen	Thymocytes and T lymphocytes
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	