

Anti-CD71 antibody [MEM-75] {PE} (STJ16100074)
STJ16100074**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD71 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-75
Concentration	
Conjugation	PE
Purification	
Dilution Range	FC using 20 µl reagent/100 µl of whole blood or 106 cells in a suspension.
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	NALM-6 human pre-B cell line
Immunogen Region	
Specificity	
Immunogen Sequence	