

Anti-PECAM-1 antibody [MEC7.46] (STJ16100083)

STJ16100083

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-PECAM-1/CD31 is suitable for use in Immunohistochemistry, ELISA, Flow Cytometry, Immunofluorescence and Immunoprecipitation research applications.
Applications	IHC-F/IHC-P/ELISA/FC/IF/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEC7.46
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:50 IF 1:50 FC 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence