

Anti-CD44 antibody [NKI-P2] (STJ16100084)

STJ16100084

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD44 is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F/IHC-P/FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NKI-P2
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:50 IF 1:50 FC 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
Uniprot ID	CD44_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	