

Anti-CD34 antibody [QBEND-10] {Biotin} (STJ16100096)

STJ16100096

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD34 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Rhesus/Cynomolgus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	QBEND-10
Concentration	0.1mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	947
Gene Symbol	CD34
Uniprot ID	CD34_HUMAN
Immunogen	Human endothelial cell membrane vesicles.
Immunogen Region	
Specificity	
Immunogen Sequence	